



Microscope Camera MMMC-A102

www.micromed.cc | info@micromed.cc

Microscope Camera MMMC-A102

Microscope Camera MMMC-A102 delivers 2.0MP resolution for detailed imaging. Built with wide-field view, capturing expansive sample areas with clarity. Features high-speed frame rates for seamless real-time observation. Integrated with advanced sensor, ensuring clear visuals with excellent color fidelity. Our HD microscope camera enhances image quality with smooth performance.

Features:

- High-definition images
- User-friendly digital interface
- Manual and auto white balance
- Lower power consumption
- Flexible zooming options

Applications:

Microscope Camera MMMC-A102 is widely employed in biological research, industrial inspection, and scientific analysis for live demonstrations, microscopy analysis and quality control.

Specifications :

Imaging Type	HD, Wide-Field
Sensor Type	CMOS Colours
Sensor Size	1/2.8?inch
Menu	All-digital UI design
Operation Mode	Mouse Operation
White Balance	Manual/Auto
Shutter Mode	Manual/Auto
Effective Pixels	2.0 MP
Frame Rate	60 fps



7136 South Yale Avenue, Tulsa OK, 74136, USA
Email: info@micromed.cc | Website: www.micromed.cc